

### THE PROFESSIONAL'S CHOICE FOR A NATURAL SMILE


Mold Chart
Anterior and Posterior Hardened Plastic Teeth




























10° Delta∆ Shallow Cusps

∆ **32X** 





32.0

Profile of 1st Premolar

8.5

Mold Articulation							
Imperial ViVO <sup>®</sup> Upper Anterior	Lower Anterior	0° Flat Plane	10° Functional	*10° <b>Delta</b> ∆ Shallow Cusps (Anatoline®)	20° Semi- Anatomical	33° Anatomical	
1G	U	33Z	34X	Δ36Χ	33XX	34M	
1L	V, VS, VL	33Z	33X, 34X	Δ34Χ, Δ36Χ	33XX	34M	
1M	W	31Z	32X, 33X	Δ32Χ, Δ34Χ	31XX	32M	
2\$	Y, YS, YL	29Z	30X	Δ30Χ	29XX	30M	
3M	W	31Z	32X, 33X	Δ32Χ, Δ34Χ	31XX	32M	
3N	Χ	31Z	32X	Δ32Χ	31XX	32M	
3\$	Y,YS, YL	29Z	30X	Δ30Χ	29XX	32M	
5L	V, VS, VL	33Z	33X, 34X	Δ34Χ, Δ36Χ	33XX	34M	
5N	Χ	31Z	32X	Δ32Χ	31XX	32M	
<b>5</b> \$	Y, YS, YL	29Z	30X	Δ30Χ	29XX	30M	
7M	W	31Z	32X, 33X	Δ32Χ, Δ34Χ	29XX	32M	
7N	Χ	31Z	32X	Δ32X	31XX	32M	
8M	W	31Z	32X, 33X	Δ32, Δ34	31XX	30M, 32M	
8N	Χ	31Z	32X	Δ32X	31XX	32M	
88	Y, YS, YL	29Z	30X	Δ30Χ	29XX	30M	
9G	U	33Z	34X	Δ36X	33XX	34M	
10M	W	31Z	32X, 33X	Δ32Χ, Δ34Χ	31XX	32M	
12M	W	31Z	32X, 33X	Δ32Χ, Δ34Χ	31XX	32M	
12N	Χ	31Z	32X	Δ32Χ	31XX	30M, 32M	
∆ <b>E12</b>	ΔΝ	31Z	32X	Δ32Χ	31XX	32M	
∆ <b>E22</b>	ΔΗ	29Z	29X	Δ30Χ	29XX	30M	
<b>∆E32</b>	ΔΗ	31Z	32X	Δ32Χ	31XX	32M	
<b>∆F42</b>	ΔΗ	30Z	32X	Δ32Χ	31XX	32M	
<b>△D55</b>	ΔΕ	30Z	29X	Δ30Χ	29XX	30M	

Bite relation will obviously alter the size requirements for lower teeth. As a general rule, the table above will provide posterior molds and sizes which are approximately correct. \*10° DeltaΔ Shallow Cusps are comparable to the Anatoline® tooth line.

#### Justi Imperial ViVO® Denture Teeth

Justi Imperial ViVO® offers traditional quality American molds using the newest Vita® shade system. The effect of this combination is outstanding esthetics and excellent functionality for the modern technician.

The ability of the teeth to blend into the adjacent dentition makes producing great looking dentures and partials much more predictable.

The cross-linking system employed in manufacturing the Imperial ViVO® tooth line further enhances resistance to wear.

## Features and Benefits of using Imperial ViVO® Denture Teeth:

American materials and manufacture:
Assuring quality and safety of materials

#### Vacuum Molded:

Denser, non-porous acrylic for longer wear

#### Hand finished:

No grinding of flash saving production time

#### **Availability:**

Through dealers or same day shipping from California factory

#### Value

Comparable in quality, esthetics and wear resistance to teeth at 3 times the price

#### **Guarantee:**

Backed by 150 years in denture tooth manufacturing. Written guarantee available upon request

#### **Mold Comparison Chart**

Imperial ViV0®	Portrait®	Imperial ViV0®	Portrait®
1 <b>G</b>	-	8M	62G
1L	12G	8N	62E
1M	13E	<b>8S</b>	76D
2\$	13D	9 <b>G</b>	21J
3M	21E	10M	45F
3N	31F	12M	65G
<b>3S</b>	42D	12N	55F
5L	11H	∆ <b>E12</b>	12E
5N	11D	<b>∆E22</b>	22E
<b>5</b> \$	22C	<b>∆E32</b>	32E
7M	45H	<b>∆F42</b>	43F
7N	75E	∆ <b>D55</b>	55D

This chart is a guide only. Please refer to actual dimensions, shape and buccal contour for optimum mold selection.

To obtain tooth order forms and instructions in the setting and processing of synthetic resin teeth please visit our website or contact:

Justi Products, Literature Dept.

A Brand Division of American Tooth Industries

1200 Stellar Dr. • Oxnard, CA 93033-2404 Toll Free (800) 235-4639 • CA (805) 487-9868 www.americantooth.com

NOTE: For achieving optimum results with our teeth, please refer to the written instructions.

#### **Shade Availability**

The Imperial ViVO® line is available in 16 of the Classical Vita® type shades including 3 bleached shades: OM1, OM2, OM3, A1, A2, A3, A3.5, A4, B1, B2, B3, B4\*\*, C1, C2, C3\*\*, C4, D2, D3. D4\*\*

\*\*These Imperial ViVO® shades can be manufactured upon request. Please inquire for quantity and terms.

#### Made in the USA

All teeth are made in the USA and backed by over one hundred fifty years of experience. All raw materials purchased from US factories.



# Outstanding esthetics and excellent functionality



#### MANUFACTURED IN THE USA BY

#### **AMERICAN TOOTH INDUSTRIES**

1200 Stellar Drive, Oxnard, CA 93033-2404 (805) 487-9868 • FAX (805) 483-8482

E-mail: info@americantooth.com · Website: www.americantooth.com

The cross-linking system employed is an interpenetrated multilayered, intersectional polymer network. 100% bondability with methyl methacrylate denture base material.

Justi® Hardened Plastic Teeth comply with ANSI/ADA Specification #15 Published in 2008, Identical adoption of ISO Standard #22112:2005

Front cover denture created by: Thomas R. Zaleske

©2016 by American Tooth Industries. All rights reserved. No part of this mold chart may be reproduced without written permission from American Tooth Industries.









