






QUALITY ORTHODONTIC SUPPLIES

VOL 6 #1













OFFERS EXPIRE 12/31/24
NOT VALID WITH ANY OTHER OFFERS
Visit: www.pozziorthodontics.com
Buy online: americantooth.com/store

CALL TODAY TO
PLACE YOUR ORDER
805-487-9868
1200 Stellar Drive
Oxnard, California 93033-2404

Pozzi Orthodontics Expanders

POZZI® HYRAX EXPANDER

The Pozzi Hyrax Expander has been designed for the US market. Made entirely of biomedical stainless steel, completely smooth and polished with no rough areas for maxium hygiene. Expander limit and directional arrow are laser etched into the surface of the expander body. Supplied with a BP Palatal Expander Key.





Packs of 1

	CODE	Ø	\leftrightarrow		0	activation turns	PRICE
		arms	body	mm	mm	for maximum expansion limit	
1 P P BP	202-C71-09	1,5	14	9	0,8	40	\$13.18
11 BP	202-C71-11	1,5	16	11	0,8	50	\$13.18
13 BP	202-C71-13	1,5	18	13	0,8	60	\$13.18

SPECIAL: BUY 200* GET 75 FREE!

SPECIAL: BUY 70* GET 25 FREE!

SPECIAL: BUY 10* GET 3 FREE!

POZZI® HYRAX EXPANDER - NO ARM

The Pozzi® Hyrax Expander – No Arm is the next product from the popular Hyrax line. Pozzi® Hyrax Expander – No Arm has special grooves to accommodate the appliance to ensure an optimal laser welding process. Expansion capacity in mm, directional arrow are laser marked on the body of the screw. Pack of 1. Comes with a Palatal Expander Key.

This also includes a .Stl file to access digitally manufactured RPEs. Please go to https://www.pozziorthodontics.com/expanders to download the file.







	CODE		\leftrightarrow		activation turns	PRICE
			body mm	mm	for maximum expansion limit	
	202-C71-09N	9 mm	14	9	40	\$10.08
T BP	202-C71-11N	11 mm	16	11	50	\$10.08
T T BP	202-C71-13N	13mm	18	13	60	\$10.08

SPECIAL: BUY 10* GET 3 FREE!

SPECIAL: BUY 70* GET 25 FREE!

SPECIAL: BUY 200* GET 75 FREE!

Easy ordering • 805-487-9868 • pozziorthodontics.com

Pozzi Orthodontics Expanders

POZZI® MICRO-EXPANDER


The Pozzi Micro-Expander are available in 8mm and 12m mm sizes. Featuring a slim design unmatched by other rapid expanders. These unique devices have 2 extra-long arms, allowing for the possibility of anterior support and ensuring excellent stability and maximum comfort in any configuration. The expander can be placed in the standard position, with a laser-marked arrow and max-expansion length visible to the user, or in the upside-down position, indicated by laser-marked arrows, to reduce the space occupied by the expander arms if the jaw is especially narrow. The Pozzi Micro-Expander can also be used to obtain dental expansion in the lower jaw. Supplied with a BP Palatal Expander Key. Packs of 1

SPECIAL: BUY 10* GET 3 FREE!

SPECIAL: BUY 70* GET 25 FREE! SPECIAL: BUY 200* GET 75 FREE!

- It can be positioned in the standard way, i.e. with the markings of the arrow and maximum expansion visible
- It can also be positioned upside down in case of very narrow palates, to minimize the bulkiness of the arms once they are bent.





POZZI® EXPANDER FOR LOWER ARCH


The Pozzi Expander for Lower Arch for Lower Arch is a specialized design to accommodate the needs of the orthodontic specialists. These Expanders may be fabricated on the lingual surface of either of the maxillary of the mandibular arch for expansion on distalization. The overall width has been reduced resulting in an improved fit, allowing for a more comfortable appliance in the patient's mouth. Supplied with a BP Palatal Expander Key. Packs of 1



SPECIAL: BUY 4* GET 1 FREE!







Easy ordering • 805-487-9868 • pozziorthodontics.com

Offer Expires December 31, 2024 – Not valid with other offers

Pozzi Orthodontics Expanders

POZZI® RAPID TELESCOPIC EXPANDER

The Pozzi Rapid Telescopic Expander guides create constant friction during the expansion movement. This expander is ideal for very narrow palates. The expanders with 9mm, 11mm and 13mm sizes feature a compact design for ease of arm bending and reducing the overall dimensions of the device. Supplied with a BP Palatal Expander Key. Packs of 1





SPECIAL: BUY 10* GET 3 FREE! SPECIAL: BUY 4* GET 1 FREE!

CODE	Ø	\leftrightarrow	~/())))))\+	Ó	activation turns	PRICE
	arms mm	body	mm	mm	for maximum expansion limit	
202-408-102	1,5	14	9	0,8	40	\$25.80
202-408-103	1,5	16	11	0,8	50	\$25.80
202-408-104	1,5	18	13	0,8	60	\$25.80

POZZI® EXPANDER WITH ORTHOGONAL ARMS

The Pozzi Expander with Orthogonal Arms is used to obtain distal, mesial, and sagittal expansion movement in the treatment of Class II and Class maloclusions. It is an ideal solution for patients needing orthopaedic maxillary expansion. Supplied with a BP Palatal Expander Key.



- Available in two expansion capacities for every therapeutic need.
- · High stability, comfort, and safety.
- Optimal positioning using orthogonal position of bending arms, reduces the device size.
- Palatal Expander Key included, allows easy intraoral activiation by the patient.
- Micro-machined guide pins providing mechanical friction throughout treatment.
- End-stop limits the possibility of over-treatment.



SPECIAL: BUY 70* GET 25 FREE! SPECIAL: BUY 10* GET 3 FREE!





Easy ordering • 805-487-9868 • pozziorthodontics.com

Offer Expires December 31, 2024 – Not valid with other offers

Sia Orthodontics Products Dist. By Pozzi Orthodontics

Leonardo Rapid Palatal Expander

LEONARDO is a rapid palatal expander with innovative features for safer, easier, and more effective use. The telescopic guides of the expander's body remain overlapped - even at its maximum opening - to ensure rigidity, stability, and expansion precision. The arms are laser welded to maximize the resistance to torsion and compression and prevent breakage. The arm's housing is designed with no cavities or blind parts, guaranteeing maximum hygiene. The mechanical locking and anti-return system prevent unintentional disassembly. International patents protect the Leonardo rapid expander.

SPECIAL
Introductory
Prices



GROUP/ITEM # DESCRIPTION 201-V30-0812 E Leonardo Rapid Palatal Expander 8mm \$51.90 201-V30-1014 E Leonardo Rapid Palatal Expander 10mm \$51.90 201-V30-1216 E Leonardo Rapid Palatal Expander 12mm \$51.90

XP Rapid Palatal Expander

The XP line of RPE, made of medical-grade stainless steel, is designed to meet the most demanding accuracy requirements, stability, and comfort. Its compact size and the chamfered shape, make it comfortable for the patient and also suitable in cases with limited space. XP's fully mechanical friction system allows expansion control avoiding the spin-back effect. The XP series doesn't use any chemical components for friction control being completely safe for the patient.

Main Features:

HOUSING OF THE ARMS IN A SINGLE PIECE BUILT: The arms are made of soft and solid steel to allow accuracy in the bending and stability in the use.

LASER WELDING OF THE ARMS: Their U-shape and the laser welding in the single-piece housing, ensure high torsion and compression resistance, eliminating risks of accidental breakage.

FUNNEL SHAPED INSERTION HOLE: The activation hole has been designed with a generous flaring funnel, to make it easier for the user to insert the expansion key.

GROUP/ITEM #	DESCRIPTION	
201-V30-0812XL	XP Rapid Palatal Expander 8mm 1-Arm	\$37.93
201-V30-1014XL	XP Rapid Palatal Expander 10mm 1-Arm	\$37.93
201-V30-1216XL	XP Rapid Palatal Expander 12mm 1-Arm	\$37.93



XP Rapid Palatal Expander 2-Arm

GROUP/ITEM #	DESCRIPTION	
201-V30-0812 X	XP Rapid Palatal Expander 8mm	\$35.95
201-V30-1014 X	XP Rapid Palatal Expander 10mm	\$35.95
201-V30-1016 X	XP Rapid Palatal Expander 12mm	\$35.95



5


Easy ordering • 805-487-9868 • pozziorthodontics.com

Sia Orthodontics Products Dist. By Pozzi Orthodontics

D-Bar

D-DGI			
GROUP/ITEM #	DESCRIPTION		CONTENT
201-B30-DB16	D-Bar - Class II Corrector - Size 16	\$235.60	<u>'</u>
201-B30-DB17	D-Bar - Class II Corrector - Size 17	\$235.60	
201-B30-DB18	D-Bar - Class II Corrector - Size 18	\$235.60	
201-B30-DB19	D-Bar - Class II Corrector - Size 19	\$235.60	
201-B30-DB20	D-Bar - Class II Corrector - Size 20	\$235.60	
201-B30-DB21	D-Bar - Class II Corrector - Size 21	\$235.60	2 Bars 2 Molar Sockets
201-B30-DB22	D-Bar - Class II Corrector - Size 22	\$235.60	1 Ruler
201-B30-DB23	D-Bar - Class II Corrector - Size 23	\$235.60	
201-B30-DB24	D-Bar - Class II Corrector - Size 24	\$235.60	
201-B30-DB25	D-Bar - Class II Corrector - Size 25	\$235.60	
201-B30-DB26	D-Bar - Class II Corrector - Size 26	\$235.60	
201-B30-DB27	D-Bar - Class II Corrector - Size 27	\$235.60	
201-B30-DB T	D-Bar – Molar Socket	\$119.80	
201-B30-DB K	D-Bar - Set 24 pcs	\$2,355.90	24 Bars (2 per size) 24 Molar Sockets 24 Rulers 1 pair of Tweezer

SPECIAL Introductory Prices



The Advantages

Control of the molar-canine segment.

The connection system has been designed with tight tolerances to increase treatment effectiveness by minimizing unwanted interarches or cuspids movements.

Effectiveness. The position of the hook on the canine's mesial side increases the appliance's effectiveness by increasing horizontal forces compared to vertical extrusion.

✓ Easy positioning. The mesiodistal and occlusal-gingival positioning lines facilitate the positioning of the socket on the tooth.

▼ Treatment time reduction. Using a D-bar improves the molar relationship making the orthodontic treatment faster.

Comfort. Its smooth and low-profile design improves patient comfort and compliance during treatment.

Morelli Dental Products Dist. By Pozzi Orthodontics

Expansion Screw Palatal Split

Hyrax type expansion screws are indicated to correct poor maxillomandibular relations. The most common clinical situations to apply are:

- Mixed/early permanent dentition cases.
- Posterosuperior bilateral crossbite, as a result of deficiency of the maxilla apical base.
- Anterior or posterior crossbites that are not too severe.

GROUP/ITEM #	DESCRIPTION	
200-65.05.011	Palatal Splint Expander 9mm	\$20.93
200-65.05.012	Palatal Splint Expander 11mm	\$20.93
200-65.05.013	Palatal Splint Expander 13mm	\$20.93

Expansion Screw Rapid Fan 10mm

The Hyrax type rapid fast fan expansion screw is intended for the anterior expansion of the upper arch, mounted in the anterior palatal region, commonly used in the treatment of patients with craniofacial anomaly (cleft palate).

GROUP/	ITEM #	DESCRIPTION

200-65.05.014 Expansion Screw Rapid Fan 10mm \$28,90

Expansion Screw Standard 11mm

Morelli expansion screws are used in fixed or mobile appliances for expansion or disjunction of the dental arches.

200-65.05.109 Expansion Screw Standard Type 11mm \$8.60

Expansion Screw Universal 9mm

Morelli expansion screws are used in fixed or mobile appliances for expansion or disjunction of the dental arches.

GROUP/ITEM	#	DESCRIPTION
------------	---	-------------

200-65.05.103 **Expansion Screw Universal 9mm** \$8,60

Expansion Screw Fan Type 8mm

The Hyrax Type rapid fast fan expansion screw is intended for the anterior expansion of the upper arch, mounted in the anterior palatal region, commonly used in the treatment of patients with craniofacial anomaly (cleft palate).

GROUP/ITEM # DESCRIPTION

200-65.05.160 Expansion Screw Rapid Fan Type 8mm \$17.23

Expansion Screw Mini 6.5mm

Morelli expansion screws are used in fixed or mobile appliances for expansion or disjunction of the dental arches.

GROUP/ITEM # DESCRIPTION

200-65,05,106 Expansion Screw Mini 6.5mm \$8.60















2024 Pozzi Ortho Mart.indd 7 9/12/24 11:11 AM

Pozzi Stainless Steel Bands

Easy Stretch Pozzi Molar Bands

The Pozzi series of molar bands are fully proportional in all dimensions; diameter, height and thickness, to provide maximum strength and to prevent cement washout due to gaping. Pozzi bands are designed with a regular temper to allow for easier "stretch" adaptation. They are available in separate lefts and rights for better seating and a closer adaptation to the natural anatomy. Pozzi bands utilize Unitek® standard sizing, making switching easy and painless. If you are tired of paying high prices for your bands, then the Pozzi line is perfect for you!

GROUP/ITEM #

DESCRIPTION

RETAIL

031-M851-808 Pozzi Band Kit - Empty \$24.63

Purchase 270 bands and receive one box FREE



SIZE	UL 29.0	UL 29.5	UL 30.0	UL 30.5	UL 31.0	UL 31.5	UL 32.0	UL 32.5			
Qty.	1	1	1	1	1	1	1	1			
SIZE	UL 33.0	UL 33.5	UL 34.0	UL 34.5	UL 35.0	UL 35.5	UL 36.0	UL 36.5			
Qty.	2	3	4	4	4	4	4	4			
SIZE	UL 37.0	UL 37.5	UL 38.0	UL 38.5	UL 39.0	UL 39.5	UL 40.0	UL 40.5			
Qty.	3	2	2	2	1	1	1	1			
SIZE	UL 41.0	UL 41.5	UL 42.0	He	nor I	~# O	uodro	nt			
Qty.	1	1	1	Upper Left Quadrant							

GROUP/ITEM # DESCRIPTION RETAIL

031-MMBK54-UL

Upper Left Quadrant Kit of 54 Bands

\$53,25

031-31.21* Upper Left Quadrant - Single *Specify size

\$.53

Otv	1	1	1	LO	wer L	.eπ Q	uaara	ınt
SIZE	LL 41.0	LL 41.5	LL 42.0	1.0		-4-0		
Qty.	3	2	2	2	1	1	1	1
SIZE	LL 37.0	LL 37.5	LL 38.0	LL 38.5	LL 39.0	LL 39.5	LL 40.0	LL 40.5
Qty.	2	3	4	4	4	4	4	4
SIZE	LL 33.0	LL 33.5	LL 34.0	LL 34.5	LL 35.0	LL 35.5	LL 36.0	LL 36.5
Qty.	1	1	1	1	1	1	1	1
SIZE	LL 29.0	LL 29.5	LL 30.0	LL 30.5	LL 31.0	LL 31.5	LL 32.0	LL 32.5

GROUP/ITEM #

*Specify size

DESCRIPTION

RETAIL

031-MMBK54-LL Lower Left Quadrant Kit of 54 Bands 031-31.11* Lower Left Quadrant - Single

\$53,25

\$.53



Upper Right Quadrant




GROUP/ITEM #	DESCRIPTION
GROUP/ITEM #	DESCRIPTIO

RETAIL \$53.25

\$.53

031-MMBK54-UR Upper Right Quadrant Kit of 54 Bands Upper Right Quadrant - Single 031-31.20*

*Specify size UR 29.0 | UR 29.5 | UR 30.0 | UR 30.5 | UR 31.0 | UR 31.5 | UR 32.0 | UR 32.5 Qty. SIZE 2 3 4 4 4 4 4 4 4 4 4 UR 37.0 UR 37.5 UR 38.0 UR 38.5 UR 39.0 UR 39.5 UR 40.0 UR 40.5 3 2 2 UR 41.0 UR 41.5 UR 42.0



GROUP/ITEM # DESCRIPTION

\$6.71

031-MMBK54-LR Lower Right Quadrant Kit of 54 Bands \$53.25 031-31.10* Lower Right Quadrant - Single

*Specify size

\$.53

Orthodontic Band Pusher & Cement Remover

Lower Right Quadrant

Stainless steel. Time saving - two tools in one. Finely serrated band pusher on one end and an all purpose scaler on the other end.

GROUP/ITEM #

Schure Dental Scaler - Band Pusher,

Double ended 6.75" (16.5 cm)

SPECIAL \$4.00



Unitek* is not a registered trademark of American Tooth Industries.

Easy ordering • 805-487-9868 • pozziorthodontics.com

8

007-8-606-349

Pozzi Stainless Steel Crowns

This stainless steel crown line is manufactured to meet the needs of all facets of dentistry. The posterior molar crowns are designed with a shallow occlusal anatomy, allowing for more healthy tooth structure after tooth reduction. The shallow occlusal cusps permit complete compatibility with opposing teeth and reduced cuspal interference during normal mastication. The crowns are parallel on the mesial and distal walls. Unlike bell-shaped crowns, these crowns will not rock or cause separation when two or more adjoining crowns are used. We have made these crowns with a thickness that prevents occlusal wear-through that sometimes occurs with the use of thinner crowns. Unlike some of the other crowns made with thinner and less durable materials, our crowns may be welded too. These are the crowns you want for making

space maintainers or Herbst appliances, or for a semipermanent restoration. Through our own research, we have determined these crowns to be of the ideal length: short enough to eliminate excess crown trimming, and long enough to provide optimum length and contour. To make our line complete to the user, we offer first and second primary molars, first and second permanent molars and first and second bicuspids. Pozzi Style 2000 crowns are a match for Unitek* Crowns. The Pozzi* Crown numbers and sizes are the same as Unitek*. Unitek* is not a trademark of American Tooth Industries.



1st Primary Molar

· · · · · · · · · · · · · · · · · · · ·	y Molai
GROUP/ITEM #	DESCRIPTION
040-40-0022	First Primary Molar Crown Kit 72 pcs., 7 sizes
040-40-22 (+size)	Reorder pack of 5. Specify crown size and quadrant. *Example: 40-22.UL1

MOLAR CROWN KIT
SPECIAL:
\$283.76
9203.70

Regular kit price \$302.04

CROWN SPECIAL: BUY 5 PKS, PAY \$16.85 PER PK, MIX & MATCH

Reorder a pack of five crowns \$20.60



2nd Primary Molar

GROUP/ITEM #	DESCRIPTION
040-40-0023	Second Primary Molar Crown Kit 72 pcs., 7 sizes
040-40-23 (+size)	Reorder pack of 5. Specify crown size and quadrant. *Example: 40-23UL1

MOLAR CROWN KIT SPECIAL: \$283.76

Regular kit price \$302.04

CROWN SPECIAL: BUY 5 PKS, PAY \$16.85 PER PK, MIX & MATCH

Reorder a pack of



1st Permanent Molar

GROUP/ITEM #	DESCRIPTION
040-40-0024	First Permanent Molar Crown Kit 72 pcs., 8 sizes
040-40-24 (+size)	Reorder pack of 5. Specify crown size and quadrant. *EExample: 40-24UL-1

MOLAR CROWN KIT SPECIAL: \$283.76

Regular kit price \$302.04

CROWN SPECIAL: BUY 5 PKS, PAY \$16.854 PER PK, MIX & MATCH

Reorder a pack of five crowns \$20.60

UR	UL
1	1
2	2
3	3
4	4
5	5
6	6
7	7
LR	LL
1	1
2	2
3	3
4	4
5	5
6	6
7	7

JR	UL
1	1
2	2
3	3
4	4
5	5
6	6
7	7
.R	LL
1	1
2	2
3	3
4	4
5	5
6	6
7	7

R	UL
1	1
2	2
2 3 4 5	2
4	4
5	4 5
6	6
7	7
3	7 LL
1	1
2	2
2 3	2 3
4	4
4 5	5
6	6
7	7
8	8

Pozzi Stainless Steel Crowns







2nd Permanent Molar

GROUP/ITEM #	DESCRIPTION
040-40-0025	Second Permanent Molar Crown Kit - 72 pcs., 8 sizes
040-40-25(+size)*	Reorder pack of 5. Specify crown size and quadrant. *Example: 40-25UL-1

MOLAR CROWN KIT SPECIAL: \$283.76

Regular kit price \$302.04

BUY 5 PKS, PAY \$16.85 PER PK, MIX & MATCH

Reorder a pack of five crowns \$20.60

1st Bicuspid

GROUP/ITEM #	DESCRIPTION
040-40-0020	First Permanent Bicuspid Crown Kit - 72 pcs., 8 sizes (uppers) 7 sizes (lowers)
040-40-20 (+size)*	Reorder pack of 5. Specify crown size and quadrant. *Example: 40-20UL1

MOLAR CROWN KIT SPECIAL: \$283.76

Regular kit price \$302.04

CROWN SPECIAL: BUY 5 PKS, PAY \$16.85 PER PK, MIX & MATCH

Reorder a pack of five crowns \$20.60

2nd Ricuspid

E DICO	JOIG
GROUP/ITEM #	DESCRIPTION
040-40-0021	Second Permanent Bicuspid Crown Kit - 72 pcs., 8 sizes (uppers) 7 sizes (lowers)
040-40-21 (+size)*	Reorder pack of 5. Specify crown size and quadrant. *Example: 40-21UL1

MOLAR CROWN KIT

Regular kit price \$302.04

CROWN SPECIAL: BUY 5 PKS, PAY \$16.85 PER PK, MIX & MATCH

Reorder a pack of five crowns \$20.60

UR UL 1 2 2 3 3 5 5 6 6 7 LL LR 2 2 3 3 5 6 6

UR	UL
00	00
0	0
1	1
2	2
3	3
4	4
5	5
6	6
LR	LL
0	0
1	1
2	2
3	3
4	4
5	5
6	6

UR	UL
0	0
1	1
2	2
3	3
4	4
5	5
6	6
7	7
LR	LL
LR O	LL 0
0	0
0	0
0 1 2	0 1 2
0 1 2 3	0 1 2 3
0 1 2 3 4	0 1 2 3 4

Pozzi[®]

These preformed components are available to the laboratory that chooses to form a custom band or crown space maintainer. You simply solder the tube to the band/crown, **Space Maintainers** insert the maintainer rest into the tubing, and crimp it at the desired point. Your doctor can even do the crimping in the office, thereby eliminating any chance of misfitting.





Narrow Space Maintainer wire insert with rests. Pkg. of 10

GROUP/ITEM # 040-40-0006 \$44.47



Narrow Space Maintainer with drop loop. Pkg. of 10 GROUP/ITEM # 040-40-0010 \$26.61



Wide Space Maintainer wire insert with rests. Pkg. of 10 GROUP/ITEM # 040-40-0007 \$44.47



Wide Space Maintainer with drop loop. Pkg. of 10 GROUP/ITEM # 040-40-0011 \$26.61



Narrow Space Maintainer wire insert without rests. Pkg. of 10 GROUP/ITEM #



Narrow Space Maintainer Tube Pkg. of 10 GROUP/ITEM #



Wide Space Maintainer wire insert without rests. Pkg. of 10 GROUP/ITEM # 040-40-0009 \$17.17



Wide Space Maintainer Tube Pkg. of 10 GROUP/ITEM # 040-40-0013 \$47.90

SPECIAL: BUY 3 GET 1 FREE - MIX&MATCH OF EQUAL OR LESSER VALUE

Easy ordering • 805-487-9868 • pozziorthodontics.com

10

Offer Expires December 31, 2024 - Not valid with other offers

Pozzi Wire

Pozzi® Retainer Wire

(Medical Specification)

Diamond drawn: yields a perfectly round wire.

- Medium hard temper: a true 1/2 hard wire.
- · Easily formed: ideal for hand forming.
- · Non-brittle: can be bent 12-14 times before breaking to prevent remakes.



MIX&MATCH &

015-19907-2 015-19907-3 015-19907-4

Size Standard Wire

DESCRIPTION Retainer Wire Straight lengths of 14", tube .025 10 pcs. \$4.92 Retainer Wire Straight lengths of 14", tube .028 10 pcs. \$4.92 Retainer Wire Straight lengths of 14", tube .030 10 pcs. \$5.41 Retainer Wire Straight lengths of 14", tube .032 10 pcs. \$5.45 Retainer Wire Straight lengths of 14", tube .036 10 pcs. \$5.45

Pozzi® Standard Round Wire

(Medical Specification)

- High Elasticity, which permits greater manipulation of this 3/4 hard wire without danger of cracking or breaking.
- · High yield point which gives the wire a greater resistance to permanent deformation.
- Uniform mechanical properties from wire to wire and batch to batch.
- Packaged in 10 14" straight lengths per tube. When ordering please specify item number and size.

GROUP/ITEM#	DESCRIPTION
015-19908-05	Standard Round Wire Size:
015-19908-07	Standard Round Wire Size:
015-19908-08	Standard Round Wire Size
015-19908-09	Standard Round Wire Size
015-19908-10	Standard Round Wire Size
015-19908-11	Standard Round Wire Size
015-19908-12	Standard Round Wire Size

e: .012 10/pkg. \$4.70 : .014 10/pkg. \$4.70 : .015 10/pkg. \$4.70 e: .016 10/pkg. \$4.70 e: .017 10/pkg. \$4.70 : .018 10/pkg. \$4.70 : .019 10/pkg. \$4.70

015-19908-13 Standard Round Wire Size: .020 10/pkg. \$4.80 015-19908-14 Standard Round Wire Size: .021 10/pkg. \$4.80 015-19908-15 Standard Round Wire Size: .022 10/pkg. \$4.80 015-19908-16 Standard Round Wire Size: .024 10/pkg. \$4.80 Standard Round Wire Size: .025 10/pkg. \$4.80 015-19908-17 015-19908-18 Standard Round Wire Size: .026 10/pkg. \$5.61 015-19908-19 Standard Round Wire Size: .028 10/pkg. \$5.61 015-19908-20 Standard Round Wire Size: .030 10/pkg. \$5.61



OF EQUAL OR LESSER

015-19908-21 015-19908-22 015-19908-23 015-20016-1

Standard Round Wire Size: .032 10/pkg \$5.61 Standard Round Wire Size: .036 10/pkg \$5.61 Standard Round Wire Size: .040 10/pkg \$5.61 Standard Round Wire Size: .045 10/pkg \$5.61 Standard Round Wire Size: .051 10/pkg \$5.61

Pozzi® Cobalt Clasp **Crozat Wire**

(Medical Specification)

- · Chromium cobalt alloy. Specifically designed for the Crozat method and for fabrication of the OSAYS appliance.
- Heat-treatable: May be soldered in open flame.
- Available in 15-23 gauge.
- Coiled spools. A softer, more malleable wire saves time when making intricate bends.



Pozzi® **Standard Round Wire**

(Medical Specification)

- 3/4 hard wire for those who want more spring to their wire.
- · Surgical-grade stainless steel.
- Spooled wire allows for increased tensile strength and less waste.



SPECIAL: BUY 3 PIECES SAME SIZE AND GET 4th PIECE FREE

(ANY SIZE - SAME PRODUCT ONLY, OF EQUAL OR LESSER VALUE)

Gauge Size	mm.	inch	ozs.#	ft.	Group/Item#	
22/23	0.6	.024	2/3	22	015-410/6	\$51.31
21	0.7	.028	3/4	22	015-410/7	\$51.31
20	0.8	.032	1	22	015-410/8	\$51.31
19	0.9	.036	11/3	22	015-410/9	\$51.31
18	1.0	.040	3/4	11	015-410/1.0	\$56.44
17	1.1	.044	1	11	015-410/1.1	\$56.44
16	1.2	.048	11/4	11	015-410/1.2	\$56.44
15	1.3	.052	11/3	11	015-410/1.3	\$56.44

\$51.31 \$51.31 \$51.31



Gauge						
Size	mm.	inch	ozs.	# ft.	Group/Item	
# 21	0.7	.028	3	90	015-400/7	\$32.60
20	0.8	.032	3	68	015-400/8	\$32.60
19	0.9	.036	3	54	015-400/9	\$32.60
18	1.0	.040	3	44	015-400/1.0	\$32.60
17	1.1	.044	3	36	015-400/1.1	\$32.60 \$32.60
16	1.2	.048	3	30	015-400/1.2	\$32.60
15	1.3	.052	3	26	015-400/1.3	\$32.60
						402100

Easy ordering • 805-487-9868 • pozziorthodontics.com

Pozzi Wire

\$23.86

Per Pka.

Pozzi®

Pozzi Super BP Ball Clasp - Series 2000

This improved ball clasp is made of high-grade #304 stainless steel wire and accepts severe bends without fracturing. The round head is approximately twice the diameter of the wire shaft and will not irritate tissue. These clasps are made from diamond-drawn wire and have a bright, clean finish. In our opinion, other ball clasps on the market today have a lower tensile strength and an oval-shaped head that does not possess adequate retention required by retainers. If you are looking for a better ball clasp to give superior retention on your removable appliances, try the 2000 series.

Available in three sizes.

GROUP/ITEM # DESCRIPTION

Ball Clasp 2000, Size .028, Pkg. 100 pcs. \$38.78 013-19814 013-19815 Ball Clasp 2000, Size .032, Pkg. 100 pcs. \$38.78

Ball Clasp 2000, Size .036, Pkg. 100 pcs. \$38.78







Pozzilov (Medical Specification)

Pozziloy retainer wire, is our own special formulation of nickel cobalt wire, available in a blue, soft, heat-treatable (1400/1600 N/sq.mm) tensile. This chromium-cobalt base alloy is highly stainless and mechanically more resistant to breakage. Each wire is color coded blue at the tip, and is 14 inches in length. This wire is made from a vacuum-fused alloy, diamond-drawn to obtain leveled and highly resistant forces. Packaged in a ten pack, it is available in .028, .030, .032 and .036.



GROUP/ITEM # DESCRIPTION

040-40-0002 Pozzilov .028 \$28.51 ea. 040-40-0003 Pozziloy .030 040-40-0004 Pozziloy .032

\$34.18 ea. 040-40-0005 Pozziloy .036



Pozzi® Half Round Wire

(Medical Specification)

Half round, semi-hard, type-304 tempered stainless steel wire to be used for reinforcement of dentures and prosthodontic appliances. Also used to fabricate special clasps. Precision cut to 12" (30.5 cm) in length. Diamond drawn for brightness, ultrasonically cleaned for dental use. Available in six (6) different sizes. Packaged in 6 pieces of 12" straight length per tube.

015-20074 12" Half Round Wire .90 mm 6/pk \$5.17 ea. 015-20075 12" Half Round Wire 1 mm 6/pk \$5.17 ea. 015-20076 12" Half Round Wire 1.2 mm 6/pk \$5.17 ea. 12" Half Round Wire 1.25 mm 6/pk \$7.47 ea. 015-20077 015-20079 12" Half Round Wire 1.7 mm 6/pk

SPECIAL: BUY 4 **GET 1 FREE! EQUAL OR LESSER VALUE MIX & MATCH**

Pozzi® Lingual Bars

Stainless steel lingual bars in four sizes: small, medium, large, extra large. More retention due to large serrated anchorage surface.

DESCRIPTION GROUP/ITEM #

Lingual Bar, Small - 55 mm. 1/pkg. 015-20081 Lingual Bar, Medium - 65 mm. 1/pkg. Lingual Bar, Large - 75 mm. 1/pkg. 015-20082 015-20083 015-20084 Lingual Bar, X-Large - 85 mm. 1/pkg.

\$3.04 ea.




Pozzi® Wire & Sheet Metal Gauge

Light-weight stainless steel. Measures round and rectangular wire with gauge and inch sizes. The only device which helps to verify sizes and maintains quality control. The wire gauge is a template which is used to measure the gauge or thickness of sheet metal or wire. Around the edge are graduated, numbered slots which fit wire or sheet metal. The metal or wire is tried in the various slots until it just fits without being forced.

GROUP/ITEM # 014-19895

DESCRIPTION Wire Gauge \$60.78



Soldering Tweezers

Stainless steel soldering tweezers with plastic handles for protection from heat of soldering. Available in two styles: straight or curved tipped.

GROUP/ITEM # **DESCRIPTION**

007-6400-1 Soldering Tweezers Straight Tip 61/2 007-6400-2 Soldering Tweezers Curved Tip 61/4

\$17.17

Pozzi **All Purpose Scissors**

7" long, stainless steel, autoclavable 1430 C, with plastic handles.

GROUP/ITEM # 007-20012

DESCRIPTION All Purpose Scissors \$37.68



30.05

Solder & Accessories



Premium Ticonium Wire

Premium quality, round clasp wire specially suited for casting with PREMIUM 100 alloy. 12" length (30cm) each wire. Comes 5 Sticks in a tube.

GROUP/ITEM #	DESCRIPTION	
015-61622	Premium Ticonimum Wire 18-gauge	\$22.84



Dental Alloy by Mesa distributed by ATI

Cobalt based A alloy for prostheses, type 5

GROUP/ITEM #	DESCRIPTION	
110-40090025A	Magnum VIP-A Cobalt Alloy, 1 kg	\$310.00



THE HISTORY OF POZZI ORTHODONTIC SILVER SOLDER Malibu



The Pozzi family has been manufacturing silver solders since 1934. From the melting of different alloys, including 65% of the purest silver, up to the making of rough bulk ingots, the quality control has been consistent, and the resulting product is remarkable. Of course, our solder does not have any cadmium or any other impurities. The rolling of the rough ingots to obtain square and round wires takes at least fifteen steps. The lamination is done with the most precise

gauged rolling mill in the industry. The high content of silver in our orthodontic solder means that very little material is needed in the soldering process, and the possibility of oxidation is nearly eliminated. Solders with a low silver content tend to oxidize heavily due to impurities. The combination of the ingredients in our silver solder and the high quality flux make this product the very best available in the world market.

Silver Solder

Pozzi® Malibu™ Orthodontic **Laboratory Silver and Brazing Solder**

CADMIUM & NICKEL FREE

A cadmium-free alloy with free flowing and ductile characteristics. Easy to use. Manufactured in the U.S.A. for the Pozzi brand.

SPECIAL PRICE		DESCRIPTION	GROUP/ITEM #
\$94.30	\$135.25	1oz Spool	031-32.50.001
\$471.15	\$676.18	5oz Spool	031-32.50.005
\$4,715.00	\$6,762.04	50oz Spool	031-32.50.050

Silver Solder A50N

With the addition of nickel the A50N offers good strength and flow characteristics and works well on stainless steel and carbide.

GROUP/ITEM #	DESCRIPTION	
A50N	Pozzi Ortho Silver Solder .032 – 1oz	\$49.42
A50N-5BLK	Pozzi Ortho Silver Solder .032 - 5oz	\$223.97





- Silver Content
- Good strength and ductility on nickel alloys to help avoid stress corrosion cracking.
 - Good color match on stainless steel.
 - Excellent fluidity that promotes penetration for good strength and elongation.

Easy ordering • 805-487-9868 • pozziorthodontics.com

13

1145° F

iNMAN Power Component Kit

COMPONENT KIT

For the fabrication of low-cost distalization appliances, the Inman Power Component sets the standard in the application and control of distalization force.

Now any quality laboratory can offer their customers appliances that are simple to adjust, completely controllable and are comfortable for patients. The Inman IPC provides all the necessary components to make an endless variety of quality appliances. With the IPC kit, 5 distalizers or a combination of distalizers and e-arches can be fabricated. The IPC's compressed coil spring provides an easy-to-adjust force controlled by a unidirectional sliding lock. The lock can be advanced 1-mm "click" at a time for accurate, predictable activation. The concept and function of the IPC is simplistic beauty and allows adjustment by the patient or patient's guardian, making it a very efficient appliance choice. In appliance fabrication, time is money, and with the Inman IPC, you will spend less time to make a quality product.

The Inman Power Component is sold as a kit, or you can purchase all the individual parts separately.

Pozzi Dental Products

Pozzi Malibu Iuman IPC Kit Inman Power Component 1 Kit/Pk.

E-ARCH

Ideal for transverse expansion of maxillary and mandibular arches where patient compliance is a concern. A tube-like framework is located on half of the arch, which is connected to the IPC on the opposing side.

IPC DISTALIZER W/ NANCE BUTTON




This appliance provides effective distal molar movement and limits mesial movement through the incorporation of a Nance button. The IPC extends through tubes to guide the maxillary first molars during distalization.

IPC DISTALIZER



USING TRANSVERSE EXPANSION

Transverse expansion and molar distalizer in one appliance. Trombone components add support during expansion and eliminate the need to cut any wires or supports before activating the Distalizer IPC components.



- The IPC is low profile and extremely comfortable compared to other components available.
- The lock is unidirectional and won't back up
- The IPC may be used for many types of movement, from simple space regainers to advanced 3-way arch development appliances.

USAGE

The infinite number of uses of this simple, dependable component will develop over time, and we encourage you to use your imagination and knowledge to help us discover all of the IPC's potential. Inman Orthodonta Laboratories will work with you on any ideas you may have for the IPC, and we will reward you for documented case results using your designs.

APPLIANCES PER KIT

APPLIANCES PER KIT

THE FORCES

The IPC compresses coil springs with an easy-to-adjust unidirectional sliding lock. The lock can be advanced 1-mm "click" at a time for accurate, predictable activation.



ACTIVATING THE IPC APPLIANCES Activation Tool



GROUP/ITEM # DESCRIPTION

031-IPC01K The IPC Kit (contains: (10) IPC Assy. (10) Welded Tubes Combination, (5) Plain Tubes 50mm, (8) Plain Tubes 4.5 mm, (1) Spool NiTi Open Coil Spring Wire \$311.58

The IPC Quantity Package:			
	031-IPC02X20	Qty 20 - IPC Inman Power Component Assy 20 Pak	\$399.36
	Refils for the IPC	Kit:	
	031-IPC02	Qty 10 – IPC Inman Power Component Assy	\$199.68
	031-IPC03	Qty 10 – Welded tube wire combination, 4.5 mm tube with 0.040 inside diameter for welding to 60 mm 0.036 wires.	\$68.64
	031-IPC04	Qty 05 - Plain tube, 50 mm in length with inside diameter of 0.040	\$13.74
	031-IPC05	Qty 08 - Plain tube, 4.5 mm with inside diameter of 0.040	\$10.00
	031-IPC06	Qty 01 - Spool of 012/.045 NiTi open coil spring wire, approx. 15.5" length	\$24.96

Easy ordering • 805-487-9868 • pozziorthodontics.com

Offer Expires December 31, 2024 – Not valid with other offers

High-Impact Acrylic Resin & Accessories

Pozzi® Ortho Resin - Cold Cure High Impact - Premium Resin - "Lo-Flo"

Pozzi Ortho Resin is a self-curing, high-impact acrylic that was developed for orthodontic use in making strong and attractive appliances and retainers. The powder is a "Lo-Flo" formulation that will rapidly absorb the monomer and will not slump, flow or run. Plastic thickness and form can be rapidly accomplished because when Pozzi Ortho powder and liquid are sprayed on the stone cast, the acrylic will stay put and even. Pozzi Ortho Resin Kit is supplied with clear powder and pink or clear liquid. Kit contains: 16 oz. powder, 8 oz. pink or clear liquid, dropper, a 3oz. jet dispenser, and instruction sheet.



Pozzi® Ortho Color Liquids

Use of these colors results in beautiful and unique appliances. Colors are dense and bright and are manufactured by our factory in Oxnard, California.



Color Guides

A display of 13 fun colors to assist patients in choosing their color.



SPECIAL \$42.16 **SPECIAL**

GROUP/ITEM #	DESCRIPTION	
041-RP6000-00	Color Shade Guide	\$50.42

			SPECIAL -	\$7.25
	GROUP/ITEM #	DESCRIPTION		
	041-RP6031-01	1 oz. Pozzi Ortho Color Liquid "C"	Pink	\$9.23
	041-RP6122-01	1 oz. Pozzi Ortho Color Liquid "A"	Clear	\$9.23
	041-RP6123-01	1 oz. Pozzi Ortho Color Liquid "B"	Trans. Pink	\$9.23
щ	041-RP6124-01	1 oz. Pozzi Ortho Color Liquid "D"	Red	\$9.23
BOTTLE	041-RP6125-01	1 oz. Pozzi Ortho Color Liquid "E"	Purple	\$9.23
	041-RP6126-01	1 oz. Pozzi Ortho Color Liquid "F"	Magenta	\$9.23
핗	041-RP6127-01	1 oz. Pozzi Ortho Color Liquid "G"	Orange	\$9.23
OUNCE	041-RP6128-01	1 oz. Pozzi Ortho Color Liquid "H"	Fluor. Yellow	\$9.23
=	041-RP6129-01	1 oz. Pozzi Ortho Color Liquid "I"	Green	\$9.23
	041-RP6130-01	1 oz. Pozzi Ortho Color Liquid "L"	Turquoise	\$9.23
	041-RP6131-01	1 oz. Pozzi Ortho Color Liquid "M"	Blue	\$9.23
	041-RP6132-01	1 oz. Pozzi Ortho Color Liquid	Black	\$9.23
	041-RP6133-01	1 oz. Pozzi Ortho Color Liquid	Gold	\$9.23

		_		-
			SPECIAL - S	\$29.28
	GROUP/ITEM #	DESCRIPTION		
	041-RP6031-08	8 oz. Pozzi Ortho Color Liquid "C"		\$49.41
	041-RP6122-08	8 oz. Pozzi Ortho Color Liquid "A"	Clear	\$49.41
	041-RP6123-08	8 oz. Pozzi Ortho Color Liquid "B"	Trans. Pink	\$49.41
щ	041-RP6124-08	8 oz. Pozzi Ortho Color Liquid "D"	Red	\$49.41
BOTTLE	041-RP6125-08	8 oz. Pozzi Ortho Color Liquid "E"	Purple	\$49.41
	041-RP6126-08	8 oz. Pozzi Ortho Color Liquid "F"	Magenta	\$49.41
- B	041-RP6127-08	8 oz. Pozzi Ortho Color Liquid "G"	Orange	\$49.41
OUNCE	041-RP6128-08	8 oz. Pozzi Ortho Color Liquid "H"	Fluor. Yellow	\$49.41
œ	041-RP6129-08	8 oz. Pozzi Ortho Color Liquid "I"	Green	\$49.41
	041-RP6130-08	8 oz. Pozzi Ortho Color Liquid "I"	Turquoise	\$49.41
	041-RP6131-08	8 oz. Pozzi Ortho Color Liquid "M"	Blue	\$49.41
	041-RP6132-08	8 oz. Pozzi Ortho Color Liquid	Black	\$49.41
	041-RP6133-08	8 oz. Pozzi Ortho Color Liquid	Gold	\$49.41

				SPECIAL - \$	24.15
	GROUP/ITEM #	DESCRIPTION			
	041-RP6031-04	4 oz. Pozzi Ortho Color Liquid "C"		Pink	\$29.48
	041-RP6122-04	4 oz. Pozzi Ortho Color Liquid "A"		Clear	\$29.48
	041-RP6123-04	4 oz. Pozzi Ortho Color Liquid "B"		Trans. Pink	\$29.48
щ	041-RP6134-04	4 oz. Pozzi Ortho Color Liquid "D"		Red	\$29.48
BOTTLE	041-RP6135-04	4 oz. Pozzi Ortho Color Liquid "E"		Purple	\$29.48
_	041-RP6136-04	4 oz. Pozzi Ortho Color Liquid "F"		Magenta	\$29.48
OUNCE	041-RP6137-04	4 oz. Pozzi Ortho Color Liquid "G"		Orange	\$29.48
9	041-RP6138-04	4 oz. Pozzi Ortho Color Liquid "H"		Fluor. Yellow	\$29.48
4	041-RP6139-04	4 oz. Pozzi Ortho Color Liquid "I"		Green	\$29.48
	041-RP6130-04	4 oz. Pozzi Ortho Color Liquid "L"		Turquoise	\$29.48
	041-RP6141-04	4 oz. Pozzi Ortho Color Liquid "M"	,	Blue	\$29.48
	041-RP6142-04	4 oz. Pozzi Ortho Color Liquid		Black	\$29.48
	041-RP6143-04	4 oz. Pozzi Ortho Color Liquid		Gold	\$29.48
_					

			SPECIAL - \$	68.14
	GROUP/ITEM #	DESCRIPTION		
	041-RP6031-16	1 pint Pozzi Ortho Color Liquid"C"	Pink	\$98.67
	041-RP6122-16	1 pint Pozzi Ortho Color Liquid"A"	Clear	\$98.67
	041-RP6123-16	1 pint Pozzi Ortho Color Liquid"B"	Trans. Pink	\$98.67
	041-RP6134-16	1 pint Pozzi Ortho Color Liquid"D"	Red	\$98.67
BOTTLE	041-RP6135-16	1 pint Pozzi Ortho Color Liquid"E"	Purple	\$98.67
3OT	041-RP6136-16	1 pint Pozzi Ortho Color Liquid"F"	Magenta	\$98.67
PINT	041-RP6137-16	1 pint Pozzi Ortho Color Liquid"G"	Orange	\$98.67
를	041-RP6138-16	1 pint Pozzi Ortho Color Liquid"H"	Fluor. Yellow	\$98.67
•	041-RP6139-16	1 pint Pozzi Ortho Color Liquid"I"	Green	\$98.67
	041-RP6130-16	1 pint Pozzi Ortho Color Liquid"I"	Turquoise	\$98.67
	041-RP6141-16	1 pint Pozzi Ortho Color Liquid"M"	Blue	\$98.67
	041-RP6142-16	1 pint Pozzi Ortho Color Liquid	Black	\$98.67
	041-RP6143-16	1 pint Pozzi Ortho Color Liquid	Gold	\$98.67

Easy ordering • 805-487-9868 • pozziorthodontics.com

15

9/12/24 11:11 AM

High-Impact Acrylic Resin & Accessories

Pozzi[®] Spruzzo[™] Colored Acrylic Powders -Solid Colors

Pozzi solid colored powders were designed for use with our clear monomer. These colors were created for those wanting a more solid color appearance. Designed for use with the salt and pepper technique, these colors are sure to catch every teenager's eye! Each color comes in a 4 ounce easy-to-use spray bottle and is available in many popular colors.



SPECIAL: \$12.42

	GROUP/ITEM #	DESCRIPTION		
뿔	040-40-0034	Spruzzo Powder 4 oz.	Yellow	\$16.39
E	040-40-0035	Spruzzo Powder 4 oz.	White	\$16.39
<u>8</u>	040-40-0036	Spruzzo Powder 4 oz.	Orange	\$16.39
ÖZ	040-40-0037	Spruzzo Powder 4 oz.	Black	\$16.39
5	040-40-0038	Spruzzo Powder 4 oz.	Brown	\$16.39
4	040-40-0039	Spruzzo Powder 4 oz.	Blue	\$16.39
	040-40-0040	Spruzzo Powder 4 oz.	Red	\$16.39
				ODEOLAL
	040-40-0040/SG	Spruzzo Shade Guide	\$31.07	SPECIAL

Pozzi[®] Spruzzo[™] Colored Acrylic Powders - Solid and Glitter Colors

Pozzi blended colored powders were designed for use with our clear monomer. These colors were created for those wanting a more solid color appearance. Designed for use with the salt-and-pepper technique, these colors are sure to catch every teenager's eye! Each color comes in a 4-ounce easy-to-use spray bottle and is available in many popular colors.



GROUP/ITEM #	DESCRIPTION			Special
040-40-0050	Spruzzo Glitter Powder 4 oz.	Yellow Glitter	\$34.39	
040-40-0051	Spruzzo Glitter Powder 4 oz.	Green Glitter	\$34.39	\$25.95
040-40-0052	Spruzzo Glitter Powder 4 oz.	Clear Glitter	\$34.39	\$20,95
040-40-0053	Spruzzo Glitter Powder 4 oz.	Hot Pink Glitter	\$34.39	

Therm-o-Flex Transparent

For Mouth Guards and Hard Bite Splints

Advantages:

Adjustable Softness **Reduced Chair Time Everlasting Fit Excellent Bond Patient Comfort** Athletic Safety

Also available in the Pink Fibered for Hybrid Dentures and to stabilize Lower Full Dentures.

GROUP/ITEM #	DESCRIPTION	
034-TF3347	Therm-0-Flex Acrylic Transparent Powder - 1 lb.	\$68.90
034-TF3305	Therm-0-Flex Acrylic Transparent Powder - 6 lb.	\$348.00
034-TF3307	Therm-0-Flex Acrylic Transparent Powder - 25 lbs.	\$1,190.00
034-TF3306	Therm-0-Flex Acrylic Liquid – 1 Quart	\$51.90
034-TF3348	Therm-0-Flex Acrylic Liquid – 11 oz.	\$32.10
034-TF3311	Therm-0-Flex Pink Fiber 1 lb.	\$86.20
034-TF3312	Therm-0-Flex Pink Fiber 6 lb.	\$375.80



Easy ordering • 805-487-9868 • pozziorthodontics.com

16

_:\$21,⁷⁵

Plastic Teeth For All Your Lab Needs

Bambino® Teeth

SPECIAL: GET 10 MIX&MATCH \$8.37 A CARD



The Bambino tooth is an aesthetic and functional solution to the problem of the empty spaces left by deciduous teeth.

Made of a copolymer cross-linked resin, this tooth line guarantees impermeability to infiltration of fluids and bacteria, thus preventing discoloration and changing of dimensional character. Studied and developed by orthodontic specialists as a simplified alternative to orthodontic appliances, these teeth are ideal for use as a space maintainer to ensure a correct line of teeth where both articulation and aesthetics are concerned. Use of these teeth in articular rehabilitations allows for proper mastication, which can only be obtained with perfect articular balance. Available in a single mold-upper and lower anteriors and upper and lower posteriors. Mounted on dark plastic blocks in 2 shades: 59 (milk-white), 61 (yellow straw-white).

GROUP/ITEM # DESCRIPTION

084-10077*	Justi Bambino Deciduous Upper Anterior Card (1x6)	\$13.50
084-10078*	Justi Bambino Deciduous Lower Anterior Card (1x6)	\$13.50
084-10177*	Justi Bambino Deciduous Upper Posterior Card (1x4)	\$13.50
084-10178*	Justi Bambino Deciduous Lower Posterior Card (1x4)	\$13.50
084-10278*	Justi Bambino Deciduous Combination Card (1x20)	\$13.50
084-10279*	Justi Bambino Deciduous Teeth Starter Kit (20 cards)	\$54.00

Pozzi® Retainer Cases

BEST VALUE!

BUY 20, & GET 25% OFF



The large and medium retainer boxes are virtually indestructible and feature "the living hinge" design with a locking mechanism. The medium case is perfect for an upper and lower retainer combination and the large size accommodates a functional type appliance. Your patients will love the combination key hole air vent. They will never lose an Expansor again. A label is provided for the patient's name, address, and phone number.

GROUP/ITEM #	DESCRIPTION	
040-65-0058LB	Large Retainer Case Blue 1/Pk. (pictured)	\$1.26
040-65-0058LP	Large Retainer Case Purple 1/Pk.	\$1.26
040-65-0058MB	Medium Retainer Case Blue 1/Pk.	\$1.22
040-65-0058MP	Medium Retainer Case Purple 1/Pk. (pictured)	\$1.22

Justi® Kura Pressure Pot

Sturdy, absolutely safe, very practical and has a large, 6-quart capacity. Made entirely of stainless steel, and the parts are available individually. This curing pot was designed to satisfy the needs of the modern technician.

GROUP/ITEM #	DESCRIPTION	
070-46-0001	Justi Kura Pot with air hose	\$470.71
070-46-0002	Air Intake Valve for Justi Kura Pot	\$7.46
070-46-004	Safety Relief Valve for Justi Kura Pot	\$27.60
070-46-0005	Pressure Release Valve for Justi Kura Pot	\$19.64
070-46-0006	Air Pressure Gauge for Justi Kura Pot	\$44.34
070-46-0007	Small Seal for Pressure Pot Safety Lock	\$5.25
070-46-0023	Air Supply Hose	\$30.20



Easy ordering • 805-487-9868 • pozziorthodontics.com

Tiger Products dist. By Pozzi Orthodontics

Deluxe Ortho Boxes

Suitable for mouth-guards, night-guards, and snoring devices. It features vent holes and a key holder. Deluxe Ortho Box is $1\frac{1}{4}$ " Deep

GROUP/ITEM # DESCRIPTION

203-D0-12MIX 12 Deluxe Ortho Box 1 ¼" \$12.93 203-D0-25* 25 Deluxe Ortho Box 1 ¼" \$23.25

Color: Black (BK), Clear (C), Forest Green (FG), Sky Blue (SB), Orange (OR), Pink (PK), Purple (PR), Red (R), Royal Blue (RB), White (W), Yellow (Y)

Ortho Boxes

An attractive container for storing retainers, mouth guards, bleaching trays, and splints. Made of flexible plastic, with a hinge design and a secure latching mechanism. These boxes are ideal for patients with removable appliances. Ortho Box is 1" Deep

GROUP/ITEM	#	DESCRIPTION
------------	---	-------------

203-0-25MIX 25 Ortho Box 1" Mixed Colors \$20.50 203-0-25* 25 Ortho Box 1" \$20.50

Color: Beige (T), Black (BK), Clear(C), Forest Green (FG), Green (G), Sky Blue (LB), Orange (O), Peach (P), Pink (PK), Purple (PR), Red (R), Royal Blue (RB), Silver (S), Violet (V), White (W), Yellow (Y)

Deep Deluxe Ortho Boxes

Made of flexible polypropylene material that will not break and resist crushing. This crown and bridge box has a positive lock-in mechanism. Top and bottom connect with a flexible hinge. It is an attractive crown and bridge tooth box that is available in a wide variety of fun colors to work with.



203-DD-25MIX 25 Deluxe Ortho Box 1 ½" Mix Colors \$30.70 203-DD-25* 25 Deluxe Ortho Box 1 ½" \$30.70

Color: Black (BK), Clear (C), Forest Green (FG), Sky Blue (SB), Purple (PR), Red (R), Royal Blue (RB), White (W), Yellow (Y)

Aligner Boxes

NEW! Slim line aligner & retainer case offers portability and convenience in a compact plastic storage box

GROUP/ITEM #	DESCRIPTION

203-A0-12MIX 12 Aligner Boxes – Plain Assorted Colors \$9.87 203-A0-2* 200 Aligner Boxes \$100.87

Color: Black (BK), Clear (C), Forest Green (FG), Green (G), Sky Blue (SB), Purple (PR), Red (R), Royal Blue (RB), White (W)

Plaster Complete Traps

Safe, reliable, and sanitary alternatives to metal plaster traps. Prevents clogged drains and expensive plumbing repairs by using disposable trap under every sink. The sanitary system catches all plasters and damaging materials heavier than water. Prevents leaks and spills. Easy to install and to change, no tools necessary. Complete Kit (Bucket with Lid, Two Fittings & Hose Set)

GROUP	/ITEM	#	DESCRIPTION
OKOU!	/	TT .	DESCRIPTION

203-DPT-3C Plaster Complete Traps 3 ½ gallon. Height-11", Diameter-12½" \$61.40
203-DPT-5C Plaster Complete Traps 5 gallons. Height-14½", Diameter-12½" \$65.78

Refill Traps

3 1/2" or 5 Gallon Bucket with Lid and Two Attached Fittings.

GROUP/ITEM #	DESCRIPTION

203-DPT-3R Disposable Refill Trap 3 ½ gallon. Height-11", Diameter-12½" \$26.32 203-DPT-5R Disposable Refill Trap 5 gallon. Height-14½", Diameter-12½" \$28.58













Easy ordering • 805-487-9868 • pozziorthodontics.com

^{*} DO-25 Choose Color

^{*0-25} Choose color

^{*}DD-25 Choose color

^{*}A0- 2 Choose Color

Thermoplastic Disc Aligners

MemoClearPro is a state-of-the-art clear thermoplastic copolyester (PET-G) specifically designed for aligners and aesthetic dental devices.

GROUP/ITEM #	DESCRIPTION	
204-3434	Thermoplastic Disc 85 x 125mm Round	\$36.58
204-3434-100	Thermoplastic Disc 85 x 125mm Round - 100	\$170.90
204-3435	Thermoplastic Disc 1.0 x 125mm Round	\$39.77
204-3435-100	Thermoplastic Disc 1.0 x 125mm Round - 100	\$183.68



Unique features

It features resilience, memory and transparency that make it unique in its kind. Its high biocompatibility makes it suitable for clinical use (aligners, splints, retainers).

Versatile

The wide range of thicknesses and accesories allows a versatile use of the MemoClearPro: aligners for dental displacements, restraints, expansions, bruxism, occlusal increases, retainers, splints, etc



Compatible

MemoClearPro fits to most thermoforming machines on the market (vacuum/pressure molding machines) and it is easy to process and compatible with acrylic and light-curing resins.

Thick mm	ness	Round Ø 125mm (4,92 inch)	Round Ø 120mm (4,72 inch)	Square 127x127mm (5x5 inch)
0,50	0.020	TP-3431	TP-3449	TP-3441
0,625	0.025	TP-3432	TP-3450	
0,75	0.030	TP-3433	TP-3439	TP-3443
0,85	0.035	TP-3434	TP-3451	
1,0	0.040	TP-3435	TP-3440	TP-3445
1,5	0.060	TP-3437	TP-3452	TP-3447
2,0	0.080	TP-3438		
Also ava	ilable pack	Pk 20 she		the item code





Half Moon

It accurately removes the plastic on the aligner for inserting button, accessories etc...

GROUP/ITEM # DESCRIPTION





Hole Key

It makes holes on the aligner for inserting elastics for 2nd/3rd Class traction, intrusion, etc...

GROUP/ITEM # DESCRIPTION

204-001-085/SM Cold Smart Plier Hole Key \$172.90





Triangular

It perfectly contours the edges of the plastic for a super-fast and precise finish of aligners

GROUP/ITEM # DESCRIPTION

204-001-087/SM Cold Smart Plier Trianglar \$172.90





Aligner Punch

It makes holes to adhere to attachments by increasing the retention of the aligner, or to insert an elastic for dental movements

GROUP/ITEM # DESCRIPTION

204-001-088/SM Cold Smart Plier Aligner Punch \$172.90





Easy ordering • 805-487-9868 • pozziorthodontics.com 20 Offer Expires December 31, 2024 – Not valid with other offers

2024_Pozzi_Ortho_Mart.indd 20 9/12/24 11:12 AM

Smart Pliers HYBRID Collection





Activation small, medium & large

The most sophisticated tools for controlling dental movements in the vestibular-lingual direction. The tiny design of Activation pliers allows the orthodontist to operate in the tiniest areas of the aligner for an unlimited series of interventions. With Small, activation points can be created in the lower incisors area; with Medium the movements in the upper laterals and the torque can be controlled; while cuspids and centrals are managed by using Large. Remember this basic rule: before creating activation points and lines, you must neutralize the black hole area by using Flattening, Spacing, with the "waxing" and "opening" method.



GROUP/IIEM #	DESCRIPTION	
204-001-074/SM	Activation Pliers - Small	\$172.90
204-001-079/SM	Activation Pliers - Medium	\$172.90
204-001-080/SM	Activation Pliers - Large	\$172.90



Vertical for rotations

A design engineered to create activation lines in a vertical direction for optimal rotations control. The couple of forces produced by Vertical, allows greater regularity and speed of movement.

GROUP/ITEM # DESCR





Torque Small

The torque control in the front dental segments in patients wearing aligners, has always been a problem for orthodontists. Now it is possible with the help of Torque Small and Large. The special tip design, creates an elongated shaped activation point for better movement control.

GROUP/ITEM #	DESCRIPTION

 204-001-082/SM
 Torque Pliers - Small
 \$172.90

 204-001-083/SM
 Torque Pliers - Large
 \$172.90



Lateral to open and close spaces

It creates activation points for dental movements in a mesial or distal direction. Lateral has a design that fits to the proximal surface of the tooth for movements optimal control.

GROUP/ITEM #	DESCRIPTION

204-001-076/SM Lateral Pliers \$172.90





Retention to make hooks

If the aligner is not stable, you don't get the patient's collaboration, the forces are not optimally transmitted, and you will have to make a new aligner. No longer now! With Retention you can create Adams clasps directly on the aligner. Choose the value corresponding to the size of the hook you want to make, on the graduated scale, and in a few seconds of thermoregulation, your aligner will be perfectly adherent to the dental arch.



GROUP/ITEM # DESCRIPTION

204-001-073/SM Retention Pliers \$ 172.90

Easy ordering • 805-487-9868 • pozziorthodontics.com



iThermal the key of the ThermalOrth Philosophy systems

SPECIAL Introductory Price

iThermal heats the pliers in just a few seconds and creates an even temperature at the tip of the tool, essential condition for a proper thermoforming.

iThermal is essential to the orthodontist who can modify aligners with simplicity, ease and precision.

GROUP/ITEM # DESCRIPTION

204-001-301 iThermal Induction Temp Control Induction \$412.62





THERMAL PROCESS: WHEN?

iThermal is essential to obtain a permanent deformation of the plastic whilst making buttons, bump planes and in operations on aligners over 1 mm thick.

iThermal technical data sheet

Width 90 mm

Length 210 mm

Height 90 mm

Net weight 0.9 kg

Gross 1.2 kg

Connection 230 v AC **voltage** 50/60 Hz

Absorption 100 W



GROUP/ITEM # 204-001-300

DESCRIPTION Timer

\$33.35

Easy ordering • 805-487-9868 • pozziorthodontics.com

0,11

Offer Expires December 31, 2024 – Not valid with other offers

Pozzi Waxes

Base Plate Wax

Three options available:

Medium Soft: dependable wax that boils out clean, reasonably priced and polishes to a high shine.

Medium Hard: same as Medium Soft but formulated with higher melting waxes and is specifically designed for warmer climates.

All Season: designed to be used most anywhere a good quality, economy wax is needed. This wax boils out clean, polishes nicely, and is formulated to securely hold teeth.

GROUP/ ITEM #	DESCRIPTION	
066-66-013/1LB	Medium Soft size sheets 76mm x 152mm x .060 thick in 1 lb. package.	\$18.20
066-66-014/1LB	Medium Hard size sheets 76mm x 152mm x .060 thick in 1 lb. package.	\$19.34
066-66-012/1LB	All Season size sheets 76mm x 152mm x .060 thick in 1 lb. package.	\$19.76
066-66-016/5LB	Medium Soft size sheets 76mm x 152mm x .060 thick in 5 lb. package.	\$42.47
066-66-017/5LB	Medium Hard size sheets 76mm x 152mm x .060 thick in 5 lb. package.	\$42.47
066-66-015/5LB	All Season size sheets 76mm x 152mm x .060 thick in 5 lb. package.	\$74.18



Justi® Wax Solvent No. 20

An efficient, safe and pleasant scented solvent used for wax elimination and shining after waxing a prosthetic case. An excellent separator for plaster, stone and any other gypsum product. Laboratories have reported to us that they also use Justi Wax Solvent #20 in the following applications:

- 1. Added to boil tank water.
- 2. To clean counter tops.
- 3. To clean tools and equipment.
- 4. To remove wax from clothing.

GROUP/ITEM #	DESCRIPT	ION
070-01-0616	1 Pint	\$21.80
070-01-0632	1 Quart	\$28.31
070-01-0699	1 Gallon	\$71.90

BUY 2 PAY SPECIAL: \$18.53 \$24.06 \$61.12 EACH



Thermoforming Material



Thermoforming blank made of PETG/EVA, including insulation foil for improved and endured transparency, for the production of:

- occlusion splints for more comfortable wearing characteristics
- hard surface
- soft inner side
- bonds with acrylate clear-transparent Featured on page 19

GROUP/ITEM #	DESCRIPTION		SPECIAL
031-D420016	Malibu Thermo Hard/Soft 120 x 2 mm	\$66.41	\$63.25
031-D420017	Malibu Thermo Hard/Soft 120 x 3 mm	\$87.73	\$83.56

Easy ordering • 805-487-9868 • pozziorthodontics.com

Impression Materials



The Hydrocolloid Alginate

WHILE THE COMPETITION CONTINUES TO TEST...
WE MAKE PERFECT IMPRESSIONS!!

HIGH PRECISION ALGINATE FOR IMPRESSIONS WITH CHROMATIC PHASE INDICATOR

120 HOURS PLUS CASTING

Thirty years ago we were the first to introduce a chromatic alginate to the United States, and now in our European facility, we are manufacturing an alginate designed for dentures, crowns, bridges and orthodontics, with these unique features:

Dimensional Stability—more than 120 hours (after washing, rinsing and removing excess water by shaking, impression must be placed in a Ziploc* baq.)

Pleasant Scent—Patients like the pleasant vanilla scent, which has a neutral taste.

Chromatic Phases:

- Violet (mixing—approx. 35 sec.)
- Rose (packing/working-approx. 40 sec.)
- Light Blue (setting approx. 45 sec.)

The impression is complete in just about 2 minutes.

- · Final impression is light blue for easy reading
- · Lead and cadmium free
- Very low dust
- · Bio-compatible
- Meets ISO 1563 and ADA #18
- Easy removal—high tear strength
- Four-pound package includes:
 4 1 lb bags and 2 measuring cups
- · 6-year freshness date
- · Antimicrobial

GROUP/ITEM #DESCRIPTION 003-2802/4 Alginmax Box of 4 - 1 lb Bags



CHROMATIC PHASES

MIXING 35 Seconds

WORKING/PACKING
40 Seconds

SETTING 45 Seconds

\$73.31

Easy ordering • 805-487-9868 • pozziorthodontics.com

24

Offer Expires December 31, 2024 – Not valid with other offers

Impression Materials



ALGINPLUS

Extra-high precision alginate for impressions with chromatic phase indicator New enriched formula for enhanced detail reproduction.

- · Stabilized hypoallergenic formula
- · Chromatic indication of the working phases for easy and precise results, regardless of water temperature and hardness
- · Dust-free powder for safe working and cleaning ease
- · Fast wetting time to obtain smooth and plastic consistency for easy mixing
- Excellent gypsum compatibility
- Very high compressive strength and recovery from deformation, exceeding ISO requirements
- Scent: tropical fruit
- · Final color: yellow
- Shelf life: 6 years
- Snap setting

2024_Pozzi_Ortho_Mart.indd 25

- · Cadmium and lead free
- Excellent biocompatibility
- · High flow and easy mixing
- Package: 453 g (1 lb)



GROUP/ITEM #DESCRIPTION Alginplus Box of 4 - 1 lb Bags 003-A6004



\$81.87

Pozzi Super-Flexible **Mixing Bowls**

Their soft, flexible texture and sturdy construction make these bowls indispensible for mixing plaster or alginate. They easily conform to the hand and make mixing a smooth and quick process. These self cleaning plasticized polyvinyl bowls are available in four convenient sizes from 350ml to the max 2250ml size.

GROUP/ITEM #	DESCRIPTION	
017-5501-10	Flexible mixing bowl small blue	\$8.74
017-5501-12	Flexible mixing bowl regular blue	\$9.40
017-5501-14	Flexible mixing bowl large blue	\$10.58
017-5501-16	Flexible mixing bowl max blue	\$17.86



Pozzi Impression Material Accessories

Labora	tory Knives		017-420
GROUP/ITEM #	DESCRIPTION		017-421
017-420	Fahnstock Large Waxing Knife - 7"	\$16.32	017-422
017-421	Fahnstock Large Waxing Knife - 5¼"	\$15.36	
017-422	Lessmann Small Waxing Knife - 5¼"	\$14.82	017-423
017-423	Lessmann Small Waxing Knife – 7"	\$15.80	
017-424	Waxing Knife with Large Spoon – 7"	\$12.89	017 /0/
017-425	Waxing Knife with Small Spoon – 7"	\$12.05	017-424
	SPECIAL: MIX & MATCH		017-425
	BUY 3 GET 1 FREE!		Gu





Pozzi Impression Material Accessories

Nylon Spatulas for Alginates



GROUP/ITEM # DESCRIPTION

003-2803/BS Alginmax Blue Spatula \$2.09 003-2803/WS Alginmax White Spatula \$2.09

SPECIAL: \$1.20 EA.

Stainless Steel Blade **Spatulas for Alginate and Plaster**

This line of spatulas was designed by dental engineers for better manipulation of alginate and plaster. The curved tips on items 017-434, 017-435 and 017-436 are the most advanced technical design. They help to mix and become a scoop, which is a great help in picking up additional material.



	GROUP/ITEM #	DESCRIPTION	
(017-430	7¼" Spatula	\$17.94
(017-431	7" Spatula	\$17.94
(017-432	8¼" Spatula	\$18.57
(017-433	8" Flexible Spatula	\$18.57
(017-434	7¾" Flexible Spatula	\$17.94
(017-435	7¾" Rigid, Heavy Blade Spatula	\$17.94
(017-436	8¾" Rigid, Heavy Blade Spatula	\$18.57
	SPE	CIAL: BUY 3 GET 1 FREE!	

Disposable Impression Trays

These Pozzi® Impression Trays are made from medical-grade plastic and can be discarded after use, or cold sterilized for repeated use. They are available as single trays in six color-coded sizes, or in kits of 12 and 60 trays. The trays feature multiple retentive grooves and are designed for orthodontic and prosthodontic use. The unique ergonomic handle is user friendly and allows for a firm grip on the tray for easier insertion and removal without slippage.



GROUP/ITEM #

DESCRIPTION

003-501-060/K

Kit of 60 assorted Duralock Color-Coded Trays, 5 pcs each size

KIT SPECIAL \$40.17

Pozzi Standard Stainless Steel Trays

Perforated stainless steel impression trays-the most popular type. Retention is assured by holes over the whole surface of the tray. A rim on the back edge stops the escape of any material, which is most unpleasant for patients. An anatomical design with rounded comfortable rims especially for adaptation to the oral cavity with respect to frenum, jaw tuberosity and retromolar pads. The dimensions come from research in laboratories throughout the world and are able to cover any anatomical situation. This means more accurate and precise impressions.




27

003-S2801-12 Kit contains one each: U1, U2, U3, U4, U5, U6, L1, L2, \$193.53

SPECIAL: \$126.79 KIT

Stainless Steel Pliers

Our studies have proven that top quality pliers will save money and produce better results. Due to the nature of orthodontic wire bending and the use of heavy wires, it is essential that durable, high-quality instruments be used. Pozzi stainless steel pliers are the best quality instruments that money can buy! When ordering, please specify item number and size.



GROUP/ITEM # DESCRIPTION

Offer Expires December 31, 2024 – Not valid with other offers

SPECIAL SALE

019-AP1034-00 3-Jaw Round Plier (Sand Finish) \$65.00

Easy ordering • 805-487-9868 • pozziorthodontics.com

Pages 28-37

Buy 3 Pliers or Instruments Same Type **Get 1 FREE SAME TYPE**

Flat Nose Pliers

General utility pliers. Heavy duty.

GROUP/ITEM # DESCRIPTION 007-8-601-018 Flat Nose Pliers 51/2" (14 cm), serrated 007-8-601-019 Flat Nose Pliers 51/2" (14 cm), smooth

007-8-601-020 Flat Nose Pliers 51/2" (14 cm), with

groove \$93.87 007-8-601-020 Flat Nose Pliers 51/2" (14 cm), with

groove \$93.87

Small Flat Nose Pliers

General utility pliers for detail work. GROUP/ITEM # DESCRIPTION

007-8-601-006 Small Flat Nose Pliers 43/4"

007-8-601-007 Small Flat Nose Pliers 43/4"



Wide Flat Nose Pliers

General utility pliers with squared beak. Heavy duty.

GROUP/ITEM # DESCRIPTION

007-8-601-001 Wide Flat Nose Pliers 51/2" (14 cm), serrated \$92.79

007-8-601-002 Wide Flat Nose Pliers 5 1/2"

(14 cm), smooth \$93.25

Round Nose Pliers

Utility pliers specifically designed with flat-edged tips and round beaks on both sides.

GROUP/ITEM # DESCRIPTION 007-8-601-022 Round Nose Pliers 51/2

(14 cm), small tip

\$90.89



Round Nose Pliers

Same pliers as 007-8-601-022 with larger beaks.

GROUP/ITEM # DESCRIPTION 007-8-601-023 Round Nose Pliers 51/2"

(14 cm), large tip



Peeso

Band and crown stretching plus contouring pliers.

GROUP/ITEM # DESCRIPTION

007-8-601-051 Peeso 6" (15 cm) \$121.59



Triplex

Combining three of the most popular contouring pliers: Johnson, Reynolds and Waldron.

GROUP/ITEM # DESCRIPTION

007-8-601-038 Triplex 51/2" (14 cm)



Special Inox Plaster Nippers

Stainless steel nippers especially handy for accurate chipping of excess plaster from laboratory work models. Cuts plaster and stone.

GROUP/ITEM # 007-444

DESCRIPTION

Special Inox Plaster Nippers

THE QUICK DESTROYER



SPECIAL SALE \$48.30



Easy ordering • 805-487-9868 • pozziorthodontics.com

29

Offer Expires December 31, 2024 - Not valid with other offers

Reynolds

Curved shortened beaks. Used in band forming and contouring.

GROUP/ITEM # DESCRIPTION 007-8-601-033 Reynolds 5" (13 cm)







Benson

Larger beaks than the Peeso. Used for band and crown stretching plus contouring.

GROUP/ITEM # DESCRIPTION 007-8-601-052 Benson 6" (15 cm)








Gordon

Curved elongated beaks. Used in band forming and contouring.

GROUP/ITEM # 007-8-601-035

DESCRIPTION Gordon 5" (13 cm)





Precision ground tips assure a firm positive grip. The narrow beaks at the tip allow accessibility to all areas.

GROUP/ITEM #

DESCRIPTION 007-8-602-000 Weingart 6" (15 cm)



How

Precision-serrated tips assure a firm, positive grip during use.

GROUP/ITEM # DESCRIPTION

007-8-602-005 How 5³/₋" (14.5 cm), straight \$91.22 007-8-602-006 How 5¹/₂" (14.0 cm), curved \$109.66



Mathieu Hemostat Forcep Needle Holders

Positive locking ratchet allows instant opening and closing. Used by many orthodontists for twisting or tying ligature wire.

GROUP/ITEM # DESCRIPTION

007-8-605-914 Mathieu Needle Holder

51/2" (14 cm), standard

007-8-605-915 Mathieu Needle Holder 51/2" (14 cm), narrow



Broken Broaches Removing Pliers

Elongated and tapered tips for removing elastics from inaccessible areas.

007-8-601-101




Anterior Band Removing Pliers

Designed for the removal of anterior bands. Prevents danger of stretching or deforming the band. Comfortable for the patient without marring teeth.

GROUP/ITEM # DESCRIPTION 007-8-602-535 Anterior Band

Removing Pliers 5" (13 cm) \$90.59



Posterior Band Removing Pliers

Extended jaw design allows the removal of posterior bands from inaccessible areas. Plastic tips are optional.

GROUP/ITEM # DESCRIPTION 007-8-602-536 Posterior Band

Removing Pliers 51/2" (13.5 cm) \$122.16



Johnson

Crimping and contouring pliers for bands, crowns and general use.

GROUP/ITEM# 007-8-601-028

DESCRIPTION Johnson 51/2" (14 cm) \$95,78



Nance/Adams Clasp Pliers

Allows the technician to produce a nance or adams clasp quickly and efficiently.

GROUP/ITEM # DESCRIPTION

007-8-601-034 Nance/Adams Clasp Pliers \$151.33



Anterior Direct Bond Bracket Remover

Wedging the tips under the bracket brace, allows for easy and efficient removal. Firm grip assures ease in horizontal use.

GROUP/ITEM # 007-8-601-074

DESCRIPTION Anterior Direct Bond **Bracket Remove** \$114.16



Hemostat / Forcep Mosquito Tip Straight

Modules sit securely in the tip. Hemostat prevents crushing of the ligature.

GROUP/ITEM # DESCRIPTION

007-8-605-005 Mosquito Tip Straight



Hemostat / Forcep Mosquito Tip Curved

Provides easy access to posterior regions.

GROUP/ITEM # DESCRIPTION

007-8-630-8014 Hemostat / Forcep Mosquito Tip

\$137.21





Hemostat has a recessed holding hook that prevents tearing of either the ligature or elastic.

GROUP/ITEM # DESCRIPTION

007-8-605-004 Hemostat / Forcep Mosquito Hook Tip



Rose Torquing Pliers

Rose Torquing Pliers wide with key .018 x .022. Designed to place torque in a specific section of archwire (single tooth) without distorting the remainder of the wire. The key is designed to torque the section between holding beaks - the torqued section will not extend past the contact points of the pliers beaks and will engage easily without transfer of torque to adjacent teeth. Labial torque and lingual torque can be placed next to each other with ease. Works with all arch wires.

GROUP/ITEM # DESCRIPTION 007-OM-805W

31

Rose Torquing Pliers Wide with key

\$126.79



Easy ordering • 805-487-9868 • pozziorthodontics.com

Offer Expires December 31, 2024 - Not valid with other offers

Pozzi Wire Benders



Use for precise bending and making loops for light wires.

GROUP/ITEM# 007-8-601-060

DESCRIPTION Andresen 51/2" (14 cm)



Adams

Squared tips. The jaws are parallel at 1 mm opening.

GROUP/ITEM# 007-8-601-043

DESCRIPTION Adams 5" (13 cm)



Angle Tweed Rectangular Arch Forming Pliers

Arch bending and torquing pliers. Non-slip parallel holding action on rectangular wire up to .025" x .208".

GROUP/ITEM # DESCRIPTION

007-8-601-076

Angle Tweed Rectangular Arch Forming Pliers 5" (13 cm)

Wire & Clasp Bending Pliers

Elongated beaks for precise bending of light wires and clasps.

GROUP/ITEM # 007-8-601-079

DESCRIPTION Wire and Clasp Bending Pliers 51/2" (14 cm)

\$102,40



Young

Wire bending pliers with 4-step beak.

GROUP/ITEM # 007-8-601-064

DESCRIPTION Young 5" (13 cm)



Has narrow concave beak on one side with three step on the other side for use on round or edgewise wire up to .020".

GROUP/ITEM #

DESCRIPTION

007-8-601-070

Tweed Loop 51/2 (14 cm)



Aderer Three Pronged Pliers

Use for opening, closing and adjusting clasps.

GROUP/ITEM # 007-8-601-040 DESCRIPTION

Aderer Three Pronged 43/4" (12 cm)



Aderer

Same pliers as 007-8-601-040 with smaller and shorter tips.

GROUP/ITEM # DESCRIPTION 007-8-601-048 Aderer Three Pronged 41/2" (11.5 cm)



Pozzi 3-Jaw Plier

For bending of Rapid Expander Arms, for wires up to 1.5 mm.

GROUP/ITEM # DESCRIPTION

019-AP1034-00 3-Jaw Round Plier (Sand Finish) \$65.00



32

Easy ordering • 805-487-9868 • pozziorthodontics.com

Offer Expires December 31, 2024 - Not valid with other offers

Pozzi Wire Benders

Aderer

Arch bending pliers. Two types - 2 groove and 3 groove.

GROUP/ITEM # DESCRIPTION

007-8-601-012

Aderer 53/4" (14.5 cm) 3 grooves

007-8-601-024 Aderer 51/2" (14 cm) 2 grooves



Arch and spring bending pliers.

GROUP/ITEM # DESCRIPTION

007-8-601-030 Waldron 5" (13 cm)



Aderer Lingual Bar Heavy Face Bow

Three-pronged box joint pliers designed for the heavy bending of lingual bars, face bows and clasps. Heavy duty.

GROUP/ITEM # DESCRIPTION

007-8-601-041 Aderer Lingual Bar

(Heavy Face Bow)

Jarabak

For precise bending of light wires. The three grooves ensure accurate bending and closing of loops.

GROUP/ITEM # DESCRIPTION

007-8-601-063 Jarabak 51/2" (14 cm)



Bird Beak Pliers

Both round and flat beaks permit the bending of rectangular and round wire up to .030"

GROUP/ITEM # DESCRIPTION 007-8-601-077 Bird Beak 5" (13

\$111.73

Hollow Chop

For contouring round arch wires. Smooth beaks.

GROUP/ITEM # 007-8-601-069

DESCRIPTION Hollow Chop



Lite Wire

Pliers designed for forming and bending up to .020" wire. Also for making and closing loops.

GROUP/ITEM # DESCRIPTION

007-8-601-078 Lite Wire 5" (13 cm)



Optical Pliers

Precision concave beaks for clasp and wire bending. Can be used on both round and edgewise wire to bend to exacting measurements.

GROUP/ITEM # 007-8-601-071

DESCRIPTION Optical 4¾" (13 cm)



Schwarz

Light wire, small optical pliers with concave beaks.

GROUP/ITEM # DESCRIPTION

007-8-601-067 Schwarz 5" (13 cm)



Easy ordering • 805-487-9868 • pozziorthodontics.com

Wire Bending Pliers

Precision round and concave beaks for clasp and wire bending. Can be used on both round and edgewise wire up to .040.

GROUP/ITEM # DESCRIPTION 007-440

Wire Bending Pliers



33

Offer Expires December 31, 2024 - Not valid with other offers

2024_Pozzi_Ortho_Mart.indd 33 9/12/24 11:12 AM

Pozzi Wire Benders & Cutters

Clasp Bending Pliers

007-442

GROUP/ITEM # DESCRIPTION Clasp Bending Pliers \$136,34



Heavy Duty Waldsachs

General forming and bending pliers. Heavy duty.

GROUP/ITEM NO. DESCRIPTION 007-443 Cutting Pliers, Waldsachs \$139,21



Heavy Duty Wire & Plate

Scissors

GROUP/ITEM # DESCRIPTION 017-437

Heavy Duty Wire and Plate Scissors



Scissors And Cutters

Universal collar and crown scissors. For cutting bands and soft wire.

GROUP/ITEM # DESCRIPTION

007-8-604-046 Crown Scissors, Smooth Edge

43/4" (12 cm), Straight

007-8-604-047 Crown Scissors, Smooth Edge 43/4" (12 cm), Curved



Lab Cutting Pliers

Designed to cut lab wires up to 0.040".

GROUP/ITEM # DESCRIPTION 007-8-603-091H Lab Cutting Pliers



Pin & Ligature **Wire Cutter**

Special carbide steel blades assure accurate, dependable cutting edges. Perfect for cutting lock and key pins. Stainless steel. Cuts soft wire up to .010". Do not use on hard wire.

GROUP/ITEM # DESCRIPTION

013-19790

Pin and Ligature Wire Cutter



Wire Cutting Pliers

Designed to cut lab wire up to 0.040".

GROUP/ITEM # DESCRIPTION

007-8-603-092 Wire Cutting Pliers \$110.23



Hemostat/Forcep Mosquito Hook tip

Hemostat has a recessed holding hook that prevents tearing of either the ligature or elastic.

GROUP/ITEM # DESCRIPTION

007-8-605-004 Hemostat / Forcep Mosquito Hook Tip



Ligature Wire Cutter

Designed to clip metal ligatures.

GROUP/ITEM # DESCRIPTION 007-8-603-093 Ligature Wire Cutter



Heavy Duty Lab Pliers

Designed to cut lab wire up to 0.032".

GROUP/ITEM # DESCRIPTION

007-8-603-096 Heavy Duty Lab Pliers



Easy ordering • 805-487-9868 • pozziorthodontics.com

34

Offer Expires December 31, 2024 - Not valid with other offers

2024_Pozzi_Ortho_Mart.indd 34 9/12/24 11:12 AM

Pozzi Classika

We wanted to manufacture a pliers like the ones you used to be able to buy in the old days - without having to mortgage the farm!

We feel that once you use one of these pliers you will know why it is a real classic!

All of the "Classika" instruments are made of corrosion resistant stainless steel. All "Classika" cutting instruments have tool steel inserts. All "Classika" cutting instruments have diamond-honed cutting edges and have a hard chrome finish. They withstand all accepted

methods of sterilization.

All "Classika" instruments have a balanced assembly wobble-free sliding lock joint. No more shim joints to deal with. All come with extended "sure grip" handles and are quaranteed to be free from manufacturing defects. Any instrument that proves defective in workmanship and materials will be repaired or replaced at our discretion free of charge. All instruments carry a two year warranty.

Black Hard Wire Cutter

Specially designed for cutting hard wires. The cutting angle is 15 degrees and the cutting edges have been diamond-honed to allow hard wire cutting. Maximum cutting capacity is .20 diameter in hard wire.

GROUP/ITEM # DESCRIPTION 040-40-0015 Black hard wire cutter \$181.19



Flush Cut Distal End Cutter with Safety Hold

The particular design of this cutter allows distal end arch wires to be cut in the patient's mouth safely. It cuts wire and then "safety holds" the loose wires.

Maximum cutting capacity: rectangular wire up to .022x.028"(0,56x0,71mm); round wire up to .020" (0,51mm); twist wire up to .016" (0,41mm) in diameter.

GROUP/ITEM # DESCRIPTION 040-40-0017R Distal end cutter

Miniature Pin & **Ligature Cutter**

This standard pin cutter is non-angled and designed to cut soft wires. This pliers is very popular for chairside work. Maximum cutting capacity is .015 in soft wire.

GROUP/ITEM # DESCRIPTION 040-40-0016

Miniature pin and ligature cutter






Micro-Miniature Light Wire Cutter with Stop Screw

This pliers is designed as a miniature ligature cutter that will cut soft wire. It easily reaches those hard-to-reach areas in the mouth. It has a factory adjusted stop screw that prevents tip damage due to overscrewing the cutter during routine procedures. Maximum cutting capacity is soft wire up to .012 in diameter.

GROUP/ITEM # DESCRIPTION 040-40-0018

Micro-Miniature cutter

15° Pin & Fine Wire Cutter



040-40-0014

15 degree pin and fine wire cutter



Easy ordering • 805-487-9868 • pozziorthodontics.com

Offer Expires December 31, 2024 - Not valid with other offers

Pozzi Tube Crimping Pliers

Stainless Steel Crown Cutting Tools



Tube Crimping Pliers

Used to crimp the tube around a space maintainer wire. It can also be used for placing stops and tubing on arch wires.

GROUP/ITEM # 040-40-0033

DESCRIPTION
Tube Crimping Pliers



Stainless Steel Crown Pliers

Small Crown Crimping Pliers

This tool is handy for anterior crowns and smaller sizes in the primary molars.

GROUP/ITEM # 040-40-0031 \$246.93



Regular Crown Crimping Pliers This is the standard in the

industry for all posterior crowns.

GROUP/ITEM # 040-40-0030 **\$246.93**



Crown Contouring Pliers

For adding contour to the crown for a better fit.

GROUP/ITEM # 040-40-0032 \$246.93



Mini-Distal End Cutter

Smaller than standard. Cuts and holds in one operation.

GROUP/ITEM # 007-8-603-090H

DESCRIPTION
Mini-Distal End Cutter
\$239.68



Easy ordering • 805-487-9868 • pozziorthodontics.com

Posterior Direct Bond Bracket Remover

Direction provides easy access to confined posterior areas. Easy vertical removal.



36

2024_Pozzi_Ortho_Mart.indd 36 9/12/24 11:12 AM

Pozzi Line Additions

Boynton-style Elastomeric Pliers

Features a smaller head with narrower arms for easier access to the posterior. An ideal alternative multiple use tool for the Crown & Bridge and Ceramics department.

GROUP/ITEM # 007-8-028-125 \$165.80



Mathieu Needle Holder

Positive locking ratchet allows instant opening and closing. Used by many orthodontists for twisting or tying ligature wire. An excellent general use instrument for Crown & Bridge and Ceramic technicians worldwide.

GROUP/ITEM # 007-8-605-915O \$127.56



V-Stop Bend Pliers

Forms a 1 mm high V bend automatically. Used to shorten arch wires or to place a strategic stop in the arch wire, such as the midline, which avoids any unwarranted displacement of the archwire. Used by laboratories for wrought wire

modifications in transitional partials.

For wires up to .025".

GROUP/ITEM # 007-8-601-045 \$121.36



Prolong Flat Nose Pliers

Serrated

Similar to a standard Mathieu except the tips are ultra fine to facilitate application of elastomeric ligature ties.

GROUP/ITEM # 007-8-605-125 \$129.52



Side Cutter Pliers - Heavy Duty

Heavy duty pliers designed to cut wire up to .040"

GROUP/ITEM # 007-8-603-096H \$221.96



Mathieu Needle Holder

Short Beak

General detail serrated utility pliers for detail work 6" - 15cm

GROUP/ITEM # 007-8-601-015 \$88.60



Side Cutter Pliers

Will cut up to .020" round wire. Features dual reciprocal springs for a lighter touch and less hand fatigue.

GROUP/ITEM # 007-8-603-098H \$112.55



Pozzi Special Guarantee

Pozzi pliers and instruments are guaranteed for a period of two years from date of purchase to be free from defects in workmanship and/or materials. Any product that proves defective in workmanship and materials will be repaired or replaced, at our discretion, free of charge.



Published by: American Tooth Industries

SEE PAGE 2-7

1200 Stellar Drive • Oxnard, CA 93033-2404 USA • WWW.pozziorthodontics.com

PRSRT STD MAIL U.S. POSTAGE

PAID

OXNARD, CA PERMIT NO. 539

NEW PRODUCTS Line of Expanders







ORTHO MATERIALS





ACCEPTED!













1200 Stellar Drive • Oxnard, California 93033-2404 | Email: info@americantooth.com | Phone 805-487-9868 | Fax (805) 483-8482